

BAUM POCKET

Aneroid sphygmomanometer



The Pocket aneroid comes from a company that has specialised in accurate sphygmomanometers for over 85 years. This elegant and compact instrument is the finest of its type on the market. The mechanisms are highly-refined and thoroughly tried and tested. Supplied with well-made, calibrated cuffs, it will give dependable, clear and accurate results for a lifetime.

- Precision low-friction movement and sophisticated twin diaphragm system
- Unique high contrast scale
- Clear pointer designed to minimise parallax error
- Attractive die cast casing
- Accurate to 1% full scale \pm 3mm Hg
- Comprehensive guarantees – see back for details

BAUM POCKET

Aneroid sphygmomanometer

Pocket Aneroid (300mm Hg)

The Baum Pocket Aneroid is a 300mm Hg light portable instrument. Like all the manometers from W.A. Baum it is built to agree with our master instruments, which are traceable to the United States National Institute of Standards and Technology. They comply with British Standards and to AAMI SP-9 for accuracy and performance (1% full scale \pm 3mm Hg). Its unique high-contrast scale and red pointer increase visibility and reduce parallax error. The manometer is manufactured with a die-cast zinc manometer housing, precision low-friction movement and twin pre-seasoned phosphor-bronze diaphragm capsules.

Every Baum Aneroid instrument is guaranteed against defects in workmanship and material for a period of one year from the date of purchase. We recommend that every aneroid instrument be checked for accuracy periodically by following our "Testing Instrument Accuracy" procedure to determine if it needs adjustment. Adjustments will be made to Baum aneroid instruments without charge when returned, postage prepaid unless the instrument has been damaged in shipment, tampered with or abused. Our pledge is that our customers will be completely satisfied with our products.

Ordering information

All of our products are supplied non latex.

For product codes for any other cuff size or special colours please see our website or call us.

Pocket Aneroid

1050-00 Pocket Aneroid with Adult Calibrated V-Lok® and case

Re-usable cuffs

The Baum Cuff is dependable, accurate and easy to use. The inner surface is marked to indicate clearly if the correct size cuff is being used, by showing the limb size range allowed. Genuine Velcro® is used for the fastener to maximize the product's life. The cuff has an inner bag or bladder and it is clearly marked on the inner and outer surface to ensure correct application on the brachial artery. The cuff material is a tightly woven, urethane coated, heat-set Dacron™ Polyester fabric, stain-resistant and dimensionally stable. All parts are non latex.

Key Features

Durable: Baum cuffs are double stitched at critical points for long wear and are treated with an effective anti-microbial agent. Cuffs, bulbs and tubing are guaranteed for five years.

Comfortable: Cuffs conform to the shape of the limb and have no hard stiff edges.

Colours: Cuffs are blue but special orders can be made in colours – blue, grey, green, brown, tan, scarlet. For example you can order a mixture of sizes of cuff in blue throughout.

Disposable cuffs for single patient use

As concerns over patient cross-contamination have grown, Baum have produced a hard wearing single-patient cuff from light blue clinical grade PVC available in four sizes. Except for the materials used, it contains the calibration and marking features of the re-usable cuff and is a high-quality product designed to give long service to a single patient.

A separate data sheet for cuffs is available.

Hawksley ● Medical and Laboratory Equipment

Marlborough Road, Lancing Business Park, Lancing, W. Sussex BN15 8TN, UK

T +44 (0)1903 752815 ● F +44 (0)1903 766050 ● www.hawksley.co.uk

