

# FREESLIDE

## Rolling pad



A revolutionary patient transfer pad which dramatically reduces the effort required to move patients laterally. The outer slip sheet slides with nearly no resistance against the board to allow a patient to be moved from bed to trolley and back.

- Widespread use for transfer of patients in the OT, ITU, recovery room and on the ward
- Folding design for storage
- Durable and washable construction
- Suitable for all body weights
- Convenient carrying handles
- 180cm x 66cm size

### Ordering information

05100-00      Freeslide patient rolling pad

**Hawksley** ● Medical and Laboratory Equipment

Marlborough Road, Lancing Business Park, Lancing, W. Sussex BN15 8TN, UK

T +44 (0)1903 752815 ● F +44 (0)1903 766050 ● [www.hawksley.co.uk](http://www.hawksley.co.uk)

