


GLASS SYRINGES




Our range of high-quality borosilicate syringes are ideal for accurate laboratory dispensation, featuring exact calibrations fused into the glass.


- Interchangeable precision-machined plungers, ground and superfinished
- Metal parts are brass, nickel then chrome plated
- Interference-fit tips with no glue. Glass-to-metal tightness guaranteed
- Borosilicate, 1st hydrolytic class, Jena type Durobax G20.
Resistant to cooling from 200–20°C in air and 100–4°C in water
- Resistant to chemical erosion and repeated sterilisation
- Quality-controlled to British Standard 1263

GLASS SYRINGES


Ordering information

Glass syringes

Description	Capacity	Code
		
Glass Tip, Luer, Centre Nozzle	0.25ml	GT14LO
Glass Tip, Luer, Centre Nozzle	0.5ml	GT12LO
Glass Tip, Luer, Centre Nozzle	1ml	GT1LOO
Glass Tip, Luer, Centre Nozzle	2ml	GT2LOO
Glass Tip, Luer, Centre Nozzle	5ml	GT5LOO
Glass Tip, Luer, Centre Nozzle	10ml	GT10LO
Glass Tip, Luer, Centre Nozzle	20ml	GT20LO
Glass Tip, Luer, Centre Nozzle	50ml	GT50LO
Glass Tip, Luer, Centre Nozzle	100ml	GT100L
Glass Tip, Luer, Centre Nozzle	150ml	GT150L




Metal Tip, Luer, Centre Nozzle	1ml	MT1LOO
Metal Tip, Luer, Centre Nozzle	2ml	MT2LOO
Metal Tip, Luer, Centre Nozzle	5ml	MT5LOO
Metal Tip, Luer, Centre Nozzle	10ml	MT10LO
Metal Tip, Luer, Centre Nozzle	20ml	MT20LO
Metal Tip, Luer, Centre Nozzle	50ml	MT50LO
Metal Tip, Luer, Centre Nozzle	100ml	MT100L
Metal Tip, Luer, Centre Nozzle	150ml	MT150L



Metal Tip, Luer Lock, Centre Nozzle	1ml	LL1LOO
Metal Tip, Luer Lock, Centre Nozzle	2ml	LL2LOO
Metal Tip, Luer Lock, Centre Nozzle	5ml	LL5LOO
Metal Tip, Luer Lock, Centre Nozzle	10ml	LL10LO
Metal Tip, Luer Lock, Centre Nozzle	20ml	LL20LO
Metal Tip, Luer Lock, Centre Nozzle	50ml	LL50LO
Metal Tip, Luer Lock, Centre Nozzle	100ml	LL100L
Metal Tip, Luer Lock, Centre Nozzle	150ml	LL150L
Metal Tip, Luer Lock, Centre Nozzle	200ml	LL200L

Tuberculin long series syringes

Description	Capacity	Code
		
Glass Tip, Luer	1ml	TB1GTO
Glass Tip, Luer	2ml	TB2GTO
Metal Tip, Luer	1ml	TB1MTO
Metal Tip, Luer	2ml	TB2MTO
Metal Tip, Luer Lock	1ml	TB1LLO
Metal Tip, Luer Lock	2ml	TB2LLO

Automatic pipetting/dispensing, borosilicate glass barrel syringes

Exact calibration. Luer lock fitting, rigid metal sleeve, pistol style handgrip, complete with tubing and pipetting cannule.



Description	Capacity	Code
Complete set as above	1ml	APDS1L
Complete set as above	2ml	APDS2L
Complete set as above	5ml	APDS5L
Complete set as above	10ml	APDS10
Complete set as above	20ml	APDS20
Spare graduated barrel for above	1ml	ADSB1L
Spare graduated barrel for above	2ml	ADSB2L
Spare graduated barrel for above	5ml	ADSB5L
Spare graduated barrel for above	10ml	ADSB10
Spare graduated barrel for above	20ml	ADSB20
Spare valve sets, all sizes		APD100

Hawksley • Medical and Laboratory Equipment

Marlborough Road, Lancing Business Park, Lancing, W. Sussex BN15 8TN, UK

T +44 (0)1903 752815 • F +44 (0)1903 766050 • www.hawksley.co.uk

