

# 5 AND 8 PLACE TOTALISERS



An intuitive differential cell counter that simplifies blood cell counts. When a total of 100 cells is counted, a bell rings. The percentage of each type of cell can then be read directly from the window above each key. A simple counter-clockwise rotation of the end knobs will zero all the counts.

- Each key has the name and picture of a particular cell above it, representing groups of leucocytes
- The key legend strips are easily removed to allow customisation
- 8 Place set to count Basophils, Eosinophils, Myelophils, Juven., Stab., Seg., Lymphocytes, Monocytes
- 5 Place set to count Basophils, Eosinophils, Seg., Lymphocytes, Monocytes
- Sturdy solid construction with firm dependable keys

## Ordering information

Z20010 8 Place Totaliser  
Z20020 5 Place Totaliser

**Hawksley** ● Medical and Laboratory Equipment

Marlborough Road, Lancing Business Park, Lancing, W. Sussex BN15 8TN, UK  
T +44 (0)1903 752815 ● F +44 (0)1903 766050 ● [www.hawksley.co.uk](http://www.hawksley.co.uk)

